

# IVOCLAR VIVADENT

## ALL CERAMIC PRODUCTS / PRODUKTE



**PRESS**  
Technology



**CAD / CAM**  
Technology



# IVOCLAR VIVADENT

## ALL CERAMIC INDICATIONS / INDIKATIONEN

	IPS Empress		IPS e.max				
	IPS Empress Esthetic	IPS Empress CAD	IPS e.max Press	IPS e.max ZirPress	IPS e.max ZirCAD	IPS e.max CAD	IPS e.max Ceram
<b>Fully anatomic / vollanatomisch</b>							
Inlays / Onlays	✓	✓					
Veneers	✓	✓		✓			
Anterior crowns / Frontzahnkronen	✓	✓					
Posterior crowns / Seitenzahnkronen	✓	✓					
<b>Layered / geschichtet</b>							
Veneers	✓	✓		✓			✓
Anterior crowns / Frontzahnkronen	✓	✓	✓		✓	✓	
Posterior crowns / Seitenzahnkronen			✓		✓	✓	
3-unit anterior bridges / 3-gliedrige Frontzahnbrücken			✓		✓		
3-unit premolar bridges / 3-gliedrige Prämolarenbrücken			✓		✓		
3-unit posterior briges / 3-gliedrige Seitenzahnbrücken					✓		
4-unit anterior briges / 4-gliedrige Frontzahnbrücken					✓		
4-unit posterior bridges / 4-gliedrige Seitenzahnbrücken					✓		
Inlay bridges / Inlaybrücken				(✓*)	✓		

\* in combination with IPS e.max ZirCAD / in Kombination mit IPS e.max ZirCAD