



**New  
Products**

# Printed Digital Models

Pow Laboratories is pleased to announce the addition of printed digital models for all Orthodontic appliances. The traditional impression taking procedure has been joined by an intra-oral partner that produces a scanned digital file. Using an STL format we can receive files from any intra-oral scanner, some examples include iTero, Trios, Sirona, Lythos or any other intra-oral scanner and we can then design from that file.

Whether you send us a traditional impression or a digital file, we have the ability to digitise your impression and import it into our software where we design either a study model or working model. We will send you the study model via e-mail for you to view and share with your patient for case analyses or send you a printed set of models. All digital files are stored at our lab. Should you request a copy or need a model printed we can do so from these files, without sending a new impression or scanned file. We are also able to digitise past study models to declutter your office space of models. Call us for more info.

Whether fixed or a removable ortho appliance, we print the models and proceed to fabricate the appliance as per your Rx.

**How to send a digital file:** either directly from your scanner selecting us from your drop down menu or if your scanner does not have the ability to send it to us directly, e-mail it to us via [shades@powlab.com](mailto:shades@powlab.com)



**Return times**  
Emailed Digital File - 1 day  
Printed Study Models- 2 days

