

Pow Laboratories Inc.

Impak Splints

Tomorrow's
Technology,
Today.



ISO 9001:2008 Certified | 2015

I Individualize with ***SPARKLES*** **GLOW IN THE DARK** **NEON** **MULTICOLOURED** Stickers and Decals to enhance patient compliance

M **Material** is a heat cured acrylic that is Strong, Comfortable, Flexible, Retentive and Bisphenol A and Pthalate free.

P **Protects** patient's teeth against Bruxing and Clenching and works to eliminate or reduce the fractures of restorations.

A **Appliance** has a **SIX Month** Warranty from date of manufacturing.

K **Knowledgeable** technical support and highly skilled staff to assist in creating an appliance that is customized and designed to meet the needs of your patients.


P **Patient Care** is equally as important as the appliance itself. All splints should be warmed under warm water for a minimum of 30 seconds prior to insertion to ensure proper fit and retention. Protective cases are provided with each appliance for proper storage when not in use. Avoid using commercial cleaners. We recommend nonabrasive toothpaste to clean splints. Finger pressure should be used to place splint in mouth. Do not have the patient bite it into position.

F **FREE** daily delivery by a Customer Relations Representative or FREE Courier Service for clients outside our delivery area.



Agnieszka Kudzia
Team Leader - Orthodontics

Fabricated in 3 days



Hassan Solker
RDT, Laboratory Manager

